

Product Datasheet

IL1RA antibody (Biotin) (orb431351)

Description	Goat polyclonal antibody against IL-1RA conjugated to Biotin
Species/Host	Goat
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant bovine interleukin-1RA
Target	IL-1RA
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.