

Product Datasheet

LEUKEMIA INHIBITORY FACTOR antibody (Biotin) (orb431357)

Description	Chicken polyclonal antibody against LEUKEMIA INHIBITORY FACTOR conjugated to Biotin
Species/Host	Gallus
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant bovine leukemia inhibitory factor
Target	LEUKEMIA INHIBITORY FACTOR
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.