

Product Datasheet

TNF ALPHA antibody (Biotin) (orb431361)

Description	Rabbit polyclonal antibody against TNF ALPHA conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant bovine tumour necrosis factor alpha (TNF- α)
Target	TNF ALPHA
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.