

Product Datasheet

IL-13 antibody (orb431382)

Description	Goat polyclonal antibody to IL-13
Species/Host	Goat
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant pig interleukin-13
Target	IL-13
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.9 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.