

Product Datasheet

IL22 antibody (Biotin) (orb431387)

Description	Rabbit polyclonal antibody against IL-22 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant pig interleukin-22
Target	IL-22
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.