

Product Datasheet

iNOS antibody (orb431392)

Description	Rabbit polyclonal antibody to iNOS
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide corresponding to aa 1131-1144 of mouse iNOS
Target	iNOS
Preservatives	0.1% Sodium Azide (NaN ₃)50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.