

Product Datasheet

GATA3 antibody (orb431396)

Description	Rabbit polyclonal antibody to GATA3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Purified partial recombinant GATA3 protein expressed in E. coli
Target	GATA3
Preservatives	0.1% Sodium Azide (NaN ₃)1 mg/ml Bovine Serum Albumin30% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.