

Product Datasheet

p27/Kip1 antibody (orb431399)

Description	Rabbit polyclonal antibody to p27/Kip1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide mapping to the carboxy terminus of human p27/Kip1
Target	p27/Kip1
Preservatives	0.1% Sodium Azide (NaN ₃)0.2% Gelatin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.