

Product Datasheet

C-JUN antibody (orb431402)

Description	Rabbit polyclonal antibody to C-JUN
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide mapping to the amino terminus of mouse c-Jun
Target	C-JUN
Preservatives	0.1% Sodium Azide (NaN ₃)0.2% Gelatin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.