

Product Datasheet

ACTIN BETA antibody (orb431406)

Description	Rabbit polyclonal antibody to ACTIN BETA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide mapping to the terminal domain of human actin beta
Target	ACTIN BETA
Preservatives	0.1% Sodium Azide (NaN ₃)0.2% Gelatin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.