

Product Datasheet

CATENIN BETA antibody (orb431407)

Description	Rabbit polyclonal antibody to CATENIN BETA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide mapping to the carboxy terminal of human catenin beta
Target	CATENIN BETA
Preservatives	0.1% Sodium Azide (NaN ₃)0.2% Gelatin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.