

Product Datasheet

CD126/IL-6R antibody (orb431438)

Description	Rabbit polyclonal antibody to CD126/IL-6R
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human CD126/IL-6R
Target	CD126/IL-6R
Preservatives	0.05%ProClin30050% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.