

Product Datasheet

C-IAP1 antibody (orb431445)

Description	Rabbit polyclonal antibody to C-IAP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human C-IAP1
Target	C-IAP1
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.