

## Product Datasheet

### GSTP1 antibody (orb431460)

|                            |                                                              |
|----------------------------|--------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GSTP1                          |
| <b>Species/Host</b>        | Rabbit                                                       |
| <b>Reactivity</b>          | Human                                                        |
| <b>Conjugation</b>         | Unconjugated                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                    |
| <b>Immunogen</b>           | Recombinant human GSTP1                                      |
| <b>Target</b>              | GSTP1                                                        |
| <b>Preservatives</b>       | 0.02% Sodium Azide50% Glycerol.<br>Phosphate buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                        |
| <b>Storage</b>             | -20°C only                                                   |
| <b>Note</b>                | For research use only                                        |
| <b>Isotype</b>             | Polyclonal IgG                                               |
| <b>Clonality</b>           | Polyclonal                                                   |
| <b>Purity</b>              | Purified                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                              |