

Product Datasheet

NEDD4 antibody (orb431488)

Description	Rabbit polyclonal antibody to NEDD4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human NEDD4
Target	NEDD4
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.