

Product Datasheet

PHOSPHOLIPASE C-GAMMA-2 antibody (orb431497)

Description	Rabbit polyclonal antibody to PHOSPHOLIPASE C-GAMMA-2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human phospholipase C-gamma-2
Target	PHOSPHOLIPASE C-GAMMA-2
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.