

Product Datasheet

AKT1 (phospho-Thr308) antibody (orb431530)

Description	Rabbit polyclonal antibody to AKT1 (phospho-Thr308)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding threonine 308 of human AKT1
Target	AKT1 (pThr308)
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.