

Product Datasheet

AKT1 (phospho-Tyr315)/AKT2 (phospho-Tyr316)/AKT3 (phospho-Tyr312) antibody (orb431532)

Description	Rabbit polyclonal antibody to AKT1 (phospho-Tyr315)/AKT2 (phospho-Tyr316)/AKT3 (phospho-Tyr312)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 315 of human AKT1/ tyrosine 316 of human AKT2/ tyrosine 312 of human AKT3
Target	AKT1 (pTyr315)/AKT2 (pTyr316)/AKT3 (pTyr312)
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.