

Product Datasheet

Bcl-xL (phospho-Ser62) antibody (orb431543)

Description	Rabbit polyclonal antibody to Bcl-xL (phospho-Ser62)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 62 of human Bcl-xL
Target	Bcl-xL (pSer62)
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.