

Product Datasheet

BRCA1 (phospho-Ser988) antibody (orb431546)

Description	Rabbit polyclonal antibody to BRCA1 (phospho-Ser988)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding serine 988 of human BRCA1
Target	BRCA1 (pSer988)
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.