

Product Datasheet

CDK6 (phospho-Tyr24) antibody (orb431557)

Description	Rabbit polyclonal antibody to CDK6 (phospho-Tyr24)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 24 of human Cdk6
Target	CDK6 (pTyr24)
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.