

Product Datasheet

ErbB3 (phospho-Tyr1328) antibody (orb431586)

Description	Rabbit polyclonal antibody to ErbB3 (phospho-Tyr1328)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 1328 of human ErbB3
Target	ErbB3 (pTyr1328)
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.