

Product Datasheet

HISTONE DEACETYLASE 2 (phospho-Ser394) antibody [Out of stock] (orb431613)

Description	Rabbit polyclonal antibody to HISTONE DEACETYLASE 2 (phospho-Ser394)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	HISTONE DEACETYLASE 2 (pSer394)
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.