

Product Datasheet

p27/Kip1 (phospho-Thr187) antibody (orb431644)

Description	Rabbit polyclonal antibody to p27/Kip1 (phospho-Thr187)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding threonine 187 of human p27/Kip1
Target	p27/Kip1 (pThr187)
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.