

Product Datasheet

VAV1 (phospho-Tyr174) antibody (orb431670)

Description	Rabbit polyclonal antibody to VAV1 (phospho-Tyr174)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 174 of human VAV1
Target	VAV1 (pTyr174)
Preservatives	0.02% Sodium Azide50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.