

Product Datasheet

ZAP 70 (phospho-Tyr493) antibody (orb431671)

Description	Rabbit polyclonal antibody to ZAP 70 (phospho-Tyr493)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Phospho specific-peptide corresponding to residues surrounding tyrosine 493 of human ZAP 70
Target	ZAP 70 (pTyr493)
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.