

Product Datasheet

CASPASE-3 antibody (orb431674)

Description	Rabbit polyclonal antibody to CASPASE-3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human caspase-3
Target	CASPASE-3
Preservatives	0.05% Proclin300 50% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.