

Product Datasheet

BLUE FLUORESCENT PROTEIN antibody (orb431704)

Description	Rabbit polyclonal antibody to BLUE FLUORESCENT PROTEIN
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	E. coli expressed recombinant blue fluorescent protein
Target	BLUE FLUORESCENT PROTEIN
Preservatives	30% Glycerol0.5% Bovine Serum Albumin5mM EDTA0.03% Proclin 300. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.