

Product Datasheet

EXODUS-2 antibody (Biotin) (orb431722)

Description	Chicken polyclonal antibody against EXODUS-2 conjugated to Biotin
Species/Host	Gallus
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant rabbit CCL21 / Exodus-2
Target	EXODUS-2
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.