

Product Datasheet

IL-15 antibody (orb431725)

Description	Rabbit polyclonal antibody to IL-15
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IHC-P, IP, WB
Immunogen	Recombinant human IL-15.
Target	IL-15
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.