

Product Datasheet

GLUTAMATE DECARBOXYLASE 1/2 antibody (orb431732)

Description	Rabbit polyclonal antibody to GLUTAMATE DECARBOXYLASE 1/2
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Synthetic peptide sequence D-F-L-I-E-E-I-E-R-L-G-Q-D-L from the C terminus of rat GAD2
Target	GLUTAMATE DECARBOXYLASE 1/2
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	65/67kDa
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.