

Product Datasheet

DYNORPHIN B antibody (orb431739)

Description	Rabbit polyclonal antibody to DYNORPHIN B
Species/Host	Rabbit
Reactivity	Porcine
Conjugation	Unconjugated
Tested Applications	IH-Fr, RIA
Immunogen	Synthetic peptide, YGGFLRRQFKVVT, corresponding to full-length porcine dynorphin B (1-13), conjugated to thyroglobulin.
Target	DYNORPHIN B
Preservatives	0.1% Sodium Azide (NaN ₃)
Form/Appearance	Serum - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(aa1-13)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Serum
Expiration Date	12 months from date of receipt.