

Product Datasheet

IL-2 antibody (orb431741)

Description	Rabbit polyclonal antibody to IL-2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant human IL-2.
Target	IL-2
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.08 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.