

Product Datasheet

NF KAPPA B INHIBITOR ALPHA antibody (orb431748)

Description	Sheep polyclonal antibody to NF KAPPA B INHIBITOR ALPHA
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, Enzyme Assay, IP, WB
Immunogen	Human NF kappa B inhibitor alpha expressed in E. coli.
Target	NF KAPPA B INHIBITOR ALPHA
Preservatives	0.1% Sodium Azide (NaN ₃)0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.