

Product Datasheet

CD184 / CXCR4 antibody (orb431766)

Description	Rabbit polyclonal antibody to CD184 / CXCR4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH-Fr, IP, WB
Immunogen	Peptide sequence corresponding to amino acids 1-14 of human CD184 / CXCR4.
Target	CD184 / CXCR4
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.