

Product Datasheet

IRAK1 antibody (orb431768)

Description	Rabbit polyclonal antibody to IRAK1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IP, WB
Immunogen	Peptide from C terminus of human IRAK
Target	IRAK1
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.