

Product Datasheet

IL-2 antibody (orb431779)

Description	Rabbit polyclonal antibody to IL-2
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant ovine IL-2.
Target	IL-2
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Serum - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Serum
Expiration Date	12 months from date of receipt.