

Product Datasheet

TGN46 antibody (orb431794)

Description	Sheep polyclonal antibody to TGN46
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Immunogen	Recombinant human TGN46.
Target	TGN46
Preservatives	0.1% Sodium Azide (NaN ₃)25% Glycerol 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Serum, diluted - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Serum
Expiration Date	12 months from date of receipt.