

Product Datasheet

TGN46 antibody (orb431795)

Description	Sheep polyclonal antibody to TGN46
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, WB
Immunogen	Recombinant human TGN46.
Target	TGN46
Preservatives	0.1% Sodium Azide (NaN ₃)0.5% Bovine Serum Albumin25% Glycerol. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.25 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.