

Product Datasheet

TGN46 antibody (orb431796)

Description

Sheep anti Human TGN46 antibody recognizes

2 (TGOLN2),...

Species/Host

Sheep

Reactivity

Human

Conjugation

Unconjugated

Tested

ICC, IF, IHC-Fr, WB

Applications
Immunogen

Recombinant human TGN46.

Target

TGN46

Preservatives

0.1% Sodium Azide (NaN₃)0.5% Bovine Serum Albumin25% Glycerol. Phosphate buffered saline

Form/Appearance

Liquid: PBS 0.1% Sodium Azide (NaN₃) 0.5% Bovine Serum Albumin 25% Glycerol

Concentration

0.25mg/ml

Storage

This product is shipped at ambient temperature. It is recommended to aliquot and store at -20°C on receipt. When thawed, aliquot the sample as needed. Keep aliquots at 2-8°C for short term use (up to 4 weeks) and store the remaining aliquots at -20°C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified

MW

110 kDa

Dilution Range

IHC-Fr: 0.1-1ug/ml, WB: .1-1ug/ml

Expiration Date

12 months from date of receipt.