

Product Datasheet

CX3CR1 antibody (orb431814)

Description	Rabbit polyclonal antibody to CX3CR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Peptide corresponding to amino acids 175 to 189 of human CX3CR1.
Target	CX3CR1
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.