

Product Datasheet

CIKS antibody (orb431839)

Description	Rabbit polyclonal antibody to CIKS
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to amino acids 554-568 of human CIKS.
Target	CIKS
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.