

Product Datasheet

GADD34 antibody (orb431876)

Description	Goat polyclonal antibody to GADD34
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide sequence AAALDLSGRRG corresponding to the C-terminus of GADD34 (NP_055145).
Target	GADD34
Preservatives	0.02% Sodium Azide 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.