

Product Datasheet

Flt3 LIGAND antibody (orb431900)

Description	Rabbit polyclonal antibody to Flt3 LIGAND
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant human Flt3-Ligand.
Target	Flt3 LIGAND
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.