

Product Datasheet

INDOLEAMINE 2,3-DIOXYGENASE antibody (orb431943)

Description	Sheep polyclonal antibody to INDOLEAMINE 2,3-DIOXYGENASE
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human indoleamine 2,3-dioxygenase.
Target	INDOLEAMINE 2,3-DIOXYGENASE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Ig Fraction - liquid
Concentration	27.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Ig Fraction
Expiration Date	12 months from date of receipt.