

Product Datasheet

IL1 BETA antibody (Biotin) (orb431950)

Description	Rabbit polyclonal antibody against IL-1 BETA conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant bovine IL-1 beta.
Target	IL-1 BETA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.25mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.