

Product Datasheet

IL-1 BETA antibody (orb431951)

Description	Rabbit polyclonal antibody to IL-1 BETA
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant bovine IL-1 beta.
Target	IL-1 BETA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.