

Product Datasheet

CD288 antibody (orb431960)

Description	Rabbit polyclonal antibody to CD288
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide corresponding to a 16 amino acid sequence within human CD288.
Target	CD288
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.