

Product Datasheet

ERK / MAPK (phospho-Thr202/pTyr204) antibody (orb431974)

Description	Rabbit polyclonal antibody to ERK / MAPK (phospho-Thr202/pTyr204)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide corresponding to an amino acid sequence within rat ERK/MAPK which includes phosphorylated Thr202 and Tyr204.
Target	ERK / MAPK (pThr202/pTyr204)
Preservatives	0.05% Sodium Azide 0.01% Bovine Serum Albumin 50% Glycerol. 10mM HEPES pH 7.5, 150mM NaCl
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.