

Product Datasheet

DYNAMIN (phospho-Ser778) antibody (orb431979)

Description	Sheep polyclonal antibody to DYNAMIN (phospho-Ser778)
Species/Host	Sheep
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide corresponding to an amino acid sequence within dynamin which includes phosphorylated Ser778.
Target	DYNAMIN (pSer778)
Preservatives	0.09% Sodium Azide 0.01% Bovine Serum Albumin 50% Glycerol. 10mM HEPES pH 7.5
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.